



Cough Detection Performance of the CoughCheck™ Algorithm in Subjects with Subacute Cough

T. deLaubenfels¹, R. Powers¹, J. Kroh¹, Y. K. Au², A. Marinovich²
¹Strados Labs, Inc, Atlanta, GA, United States, ²Strados Labs Inc, Philadelphia, PA, United States
This work was sponsored by Strados Labs

RATIONALE

- Progress in deep neural network algorithms presents an opportunity to expand the application of objective cough monitoring by supplementing manual human annotation
- Cough annotation by algorithm from sound recordings would expand the scope and lower the cost of objective cough monitoring, particularly in studies with larger populations
- We report here the performance of the **CoughCheck™ signal processing algorithm** to identify cough events in continuous, real-world respiratory recordings in a population with subacute cough

LEXICON

- Precision**
Ratio of True Positive cough predictions over the sum of coughs identified by the algorithm
- Sensitivity**
Ratio of True Positive cough predictions over the sum of human-annotated coughs
- Absolute cough count**
Number of algorithm-predicted coughs vs. the number of human-annotated coughs, independent of timing
- Temporal detection of cough events**
The timing-dependent assessment of algorithm-predicted cough events vs. human-annotated cough events

METHODS

- CoughCheck was trained to identify cough events in a set of respiratory sound recordings
- Collected using a dedicated biosensor (**Fig. 1**)
 - Collected during ordinary daily conditions (including sleep, sedentary time, and activity) with **no specific restrictions** on environment (e.g., avoiding high-noise environments)
 - 249** subjects with cough (chronic (≥8 weeks) and subacute cough (3-8 weeks))
 - 3,099** hours (**233,626** human-annotated coughs)
 - CoughCheck agreement with cough identification with human-annotation was assessed in a separate benchmark set



Figure 1: The RESP Biosensor

RESULTS

Performance Benchmark Set

- Data Collection**
- Dedicated biosensor (**Fig. 1**)
 - Continuous & ambulatory
 - Ordinary daily conditions

N = 36

Condition
Subacute Cough
(3-8 weeks)

Average Duration
23.5 hours

Total
845 total hours
14,486 human-annotated coughs

Mean Cough frequency
17.11 ± 23.24 coughs/h (arithmetic)
8.53 ± 3.34 coughs/h (geometric)

24-hour Cough Frequency Determination

- Correlation with human annotation (**Fig. 2**): **R²=0.993**
- Mean absolute error (**Fig. 2**): **1.16** coughs/h
- Bland Altman analysis (**Fig. 3**):
 - Mean difference (bias): **-0.285** coughs/h
 - Standard deviation of differences: **1.77** coughs/h




Figure 2. Algorithm-predicted vs human-annotation derived hourly cough frequency




Figure 3. Bland Altman Plot of Algorithm vs. Human - Derived Cough Frequency

Precision & Sensitivity

- Absolute cough count**
- Precision **97.4%**
 - Sensitivity **95.7%**
 - Mean overcounting rate: **0.433** coughs/h

- Temporal cough events count (Fig. 4)**
- Precision: **91.5%**
 - Precision at 80% sensitivity: **90.4%**
 - Precision at 90% sensitivity: **78.8 %**




Figure 4: Precision/Sensitivity Curve for Cough Detection (temporal events)

CONCLUSIONS

- CoughCheck was able to derive 24-hour cough frequencies with satisfactory agreement with human annotation and a low mean overcounting rate in a population with subacute cough
- These results warrant the continuing development and validation of CoughCheck for objective cough monitoring purposes
- Future work will evaluate the performance of CoughCheck in other cough categories (e.g. chronic cough), etiologies and applications